


Reading graduated cylinders


Read the bottom center of each concave meniscus. Estimate halfway between minor marks when needed.

1. 10 mL / 0.1 mL marks




Reading: _____ mL

2. 10 mL / 0.1 mL marks




Reading: _____ mL

3. 10 mL / 0.1 mL marks




Reading: _____ mL

4. 10 mL / 0.1 mL marks




Reading: _____ mL

5. 10 mL / 0.1 mL marks




Reading: _____ mL

6. 10 mL / 0.1 mL marks




Reading: _____ mL

7. 10 mL / 0.1 mL marks




Reading: _____ mL

8. 10 mL / 0.1 mL marks




Reading: _____ mL

9. 10 mL / 0.1 mL marks




Reading: _____ mL

10. 10 mL / 0.1 mL marks




Reading: _____ mL

11. 10 mL / 0.1 mL marks



Reading: _____ mL

12. 10 mL / 0.1 mL marks



Reading: _____ mL