


Reading graduated cylinders


Read the bottom center of each concave meniscus. Estimate halfway between minor marks when needed.

1. 10 mL / 0.1 mL marks




Reading: _____ mL

2. 10 mL / 0.1 mL marks




Reading: _____ mL

3. 10 mL / 0.1 mL marks




Reading: _____ mL

4. 10 mL / 0.1 mL marks




Reading: _____ mL

5. 10 mL / 0.1 mL marks




Reading: _____ mL

6. 10 mL / 0.1 mL marks




Reading: _____ mL

7. 10 mL / 0.1 mL marks




Reading: _____ mL

8. 10 mL / 0.1 mL marks




Reading: _____ mL

9. 10 mL / 0.1 mL marks




Reading: _____ mL

10. 10 mL / 0.1 mL marks




Reading: _____ mL

11. 10 mL / 0.1 mL marks



Reading: _____ mL

12. 10 mL / 0.1 mL marks




Reading: _____ mL

Reading graduated cylinders - Answer key


Read the bottom center of each concave meniscus. Estimate halfway between minor marks when needed.

1. 10 mL / 0.1 mL marks




Reading: 7.75 mL

2. 10 mL / 0.1 mL marks




Reading: 9.8 mL

3. 10 mL / 0.1 mL marks




Reading: 1.1 mL

4. 10 mL / 0.1 mL marks




Reading: 3.95 mL

5. 10 mL / 0.1 mL marks




Reading: 3.9 mL

6. 10 mL / 0.1 mL marks




Reading: 6.35 mL

7. 10 mL / 0.1 mL marks




Reading: 8.75 mL

8. 10 mL / 0.1 mL marks




Reading: 2 mL

9. 10 mL / 0.1 mL marks




Reading: 3.55 mL

10. 10 mL / 0.1 mL marks




Reading: 6.8 mL

11. 10 mL / 0.1 mL marks



Reading: 8.95 mL

12. 10 mL / 0.1 mL marks



Reading: 9.85 mL